

Fedegari Vertical Laboratory Sterilisers



FVA Small Laboratory Sterilisers

- Compact and Ergonomic Design
- High Performance and Versatility
- Modular Design and Construction
- Optimised for Energy Saving
- Extremely Safe and Reliable
- Fully Validated Software

Pharma innovation comes to your lab



Standard build includes

- DCSPLUS10 Colour Touch Screen Process Controller
- Electro-polished AISI316L Chamber
- Exclusive Heat Resistant Maintenance Free Door Gasket
- 30 Configurable Sterilisation Cycles
- External Paneling wholly in AISI304 Stainless Steel



Commonly selected options

- Integral Electric Steam Generator
- Air Steam Mixture (Counter pressure) Functionality
- Pathogen Decontamination (up to CAT3)
- Load and Effluent Cooling
- EN285/HTM2010 Compliance
- 21CFRpt11 - Data storage, Logging and Electronic Signatures

Many more options are available. These can be incorporated during the initial build or owing to the 'fully modular' design can be retrofitted at a later date.

Facts and Figures

	FVA2/A1	FVA3/A1
Useful capacity (l)	75	140
Chamber diameter (mm)	400	500
Useful height (mm)	600	700
Ext dimensions – W	1200	1200
Ext dimensions - H	1305	1305
Ext dimensions - D	875	875
Steam gen (kW)	6.5	7.5

Contact sales@steritech.co.uk for a quotation or Telephone +44 (0)1293 824400 (ext3)

SteriTech

Stilmas Equipment

- Pre-treatment
- Purified Water
- Multi Effect Still
- Vapour Compression Still
- Pure Steam Generators
- PS & WFI Combinations
- Storage & Distribution Systems

Fedegari Equipment

- Saturated Steam Sterilisers
- Counterpressure Sterilisers
- UltraWashers
- Dry Heat Sterilisation
- Pathogen Sterilisers
- Combined Units
- Laboratory Sterilisers

Fedegari non-thermal equipment

- Low Temperature Surface Decontamination Stations
- Laminar Flow Carts
- Isolators and Cleanroom Contamination Controls

Olsa Equipment

- Solids: Filtering, mixing and drying
- Semi-Solids: Mixing and turnkey plants
- Liquids: Preparation vessels, fermenters and bio reactors

Bano Equipment

- Stainless Steel Cleanroom Furniture
- Specialist Cleanroom Equipment

SteriTech Limited

- Complete Turn-Key Projects